

# New Clamps/Valves



## ACF Big Pharma Clamp

### Part No. ACF-C100

New reusable Pharma Clamp for tubing up to 1" ID (25,4 mm)

Dim. 77x65x19 mm; Seat distance when closed 3,6 mm, tubing wall should be at least 2,38 mm (0,094").

## ACF Secure Valve Clamp

### Part No. ACF-C075

New reusable Secure Valve Clamp for tubing up to 0,700" OD (18,0 mm)

Dim. 70x30x20 mm; Will operate up to 100 psi when placed correctly on the tubing.



Material: Borosilicate (glass) reinforced nylon.

Sterilizable by Autoclave, steam, irradiation.

\*Warning!

Maximum irradiation dosage is 6 mega rads, dosage is cumulative. Over this value you could experience discolouration and brittleness.

Do not over torque, could damage tubing and valve. Maximum torque is 2,0 Nm.

Order No.	Description	For tubing Size (mm)	For tubing Size (inch)	Dimensions (mm)	Material	Colour
<b>ACF-C100</b>	Big Pharma-Clamp	25,4	1	77x65x19	Borosilicate reinf. nylon	white
<b>ACF-C075</b>	Secure Valve-Clamp	18,0	0,7	70x30x20	Borosilicate reinf. nylon	white
<b>DE1</b>	Tube Clamp	0—17	0-0,65	45x38	Delrine	Yellow/red
<b>PS1</b>	Tube Clamp	0—17	0-0,65	45x38	Polysulfone	Yellow/pink
<b>INF14007</b>	Pinch Clamp	6,8	0,268	39x26	PP	white
<b>INF11014</b>	Slide Clamp	16,0	0,375	72x32	HDPE	white
<b>INF01040</b>	Pinch Clamp	6,0	0,2362	24x15	PP	nature
<b>INF04060</b>	Pinch Clamp	8,1	0,3188	30x18	Polyacetal	nature
<b>INF02060</b>	Pinch Clamp	18,8	0,433	62x34	Polyacetal	white
<b>INF340TCL</b>	Medi Clamp Large	18,8	0,433	56x34	PE	nature